



Temporary Bathroom / Toilet - Comparison Chart

Purpose: this list provides a starting point for assessors who are looking at equipment options from a range of suppliers. Assessors should contact suppliers for full details and specifications.

This list includes some options requested through Enable New Zealand and does not include all options available in NZ. Information is current as of the date shown in the footer and subject to change.

Considerations for selection of a temporary shower and/or toilet:

- All orders should be submitted as part of the ACC Housing Modifications service (<https://www.acc.co.nz/assets/contracts/mres-og.pdf>)
- All potential options should be considered and documented on a Housing Modification Assessment report (ACC257 or ACC263). For example, natural supports, short term solutions, adaptive equipment
- Options presented in this list are not suitable as permanent solutions
- When considering options where maintenance is required, ensure the client/family have the capacity to do this. For example, refilling gas bottles, emptying toilet cassette etc.
- Both hire and purchase options should be explored for cost comparison. It may be better value for money to purchase rather than hire a temporary bathroom for medium-term use, as a purchased unit can be reissued. For example, extensive permanent modifications may take more than 18months to complete, so purchase may be more cost effective overall.

Supplier	Model	Geographical cover	Hire	Purchase	Reissue	Specifications	Practical considerations
Louna Bathrooms https://www.louna.org.nz/	Louna	All North Island	Yes	Yes	No	<ul style="list-style-type: none"> • 3mL x 1.8mW, 2.4mH (external) • Includes a 2m ramp (1:12 gradient) 	<ul style="list-style-type: none"> • Contains accessible shower, toilet & basin • Addresses shorter-term needs • Towable • Fits into smaller spaces - can be hiab/craned into position if required • Cost includes 6-weekly service of unit i.e. macerator pump, gas bottles filled • Need access to gully trap & outside tap • Power from extension lead or caravan plug outside home • Gas for water heating
Outback Bathrooms https://www.outbackbathrooms.co.nz/wheelchair-accessible-portable-bathrooms	Lowa Towa	All NZ	Yes	No	N/A	<ul style="list-style-type: none"> • Internal – 2.30mL x 1.45mW x 1.95mH • A clip-on, 2.4m-long wheelchair ramp is available 	<ul style="list-style-type: none"> • Contains accessible shower, toilet & basin • Addresses shorter-term needs • Towable • Solar powered (can be topped up from mains power if needed) • Gas hot water • Cassette toilet • Client responsible for refilling gas bottles & emptying toilet cassette • Garden hose connection
Highlander	Bathroom/sleepout	All NZ	No	Yes	Yes	<ul style="list-style-type: none"> • 8mL x 3mW • Shower 1.8m x 1.2m 	<ul style="list-style-type: none"> • Contains accessible shower, toilet & basin with accommodation/bedroom • Addresses medium-term needs

Highlander bathroom/sleepout							<ul style="list-style-type: none"> • Hiab/craned into position • Client responsible for servicing of unit, refilling gas bottles etc. • Power from extension lead or caravan plug outside home • Gas for water heating • Plumbing, all internal and under floor waste, entry and exit points
Highlander Highlander Accessible Bathroom	Bathroom	All NZ	No	Yes	Yes	<ul style="list-style-type: none"> • 3.66L x 2.5mW • Shower 1.8m x 1.2m 	<ul style="list-style-type: none"> • Contains accessible shower, toilet & basin • Addresses medium-term needs • Hiab/craned into position • Client responsible for servicing of unit, refilling gas bottles etc. • Power from extension lead or caravan plug outside home • Gas for water heating • Plumbing, all internal and under floor waste, entry and exit points
Highlander Highlander Pop-up Bathroom	Pop-up bathroom	All NZ	No	Yes	Yes	<ul style="list-style-type: none"> • 2.93L x 2.13mW x 2.3mH (with roof up) • Shower 1.8m x 1.2m 	<ul style="list-style-type: none"> • Contains accessible shower, toilet & basin • Addresses shorter-term needs • Designed to sit in an internal garage • Placed on the floor of the garage and connected to the closest power and water sources • Inbuilt macerator pump to transfer waste to the nearest toilet or holding tank
Highlander Highlander Mini Bathroom	Highlander Mini	All NZ	No	Yes	Yes	<ul style="list-style-type: none"> • 1.2x1.2m 	<ul style="list-style-type: none"> • Accessible shower only, toilet only or both • Addresses shorter-term needs • Shower dome • Two-stage shower curtain • Can be broken down into three main parts and installed anywhere e.g. bedroom/lounge • Connects to water and waste
Active Healthcare https://activehealthcare.co.nz/products/hygiene/careport/careport-freestanding-shower-cubicle-extra/	Careport Extra ACC Complex List	All NZ	No	Yes	Yes	Shower only – 1.2mW x 1.4mD x 2.1mH	<ul style="list-style-type: none"> • Freestanding portable shower unit • Toilet and/or basin can be added • Addresses shorter-term needs • Can be installed on almost any floor surface • Connected to existing plumbing or temporarily connected to an outside water source • Pumped waste options to suit drainage conditions
Active Healthcare https://activehealthcare.co.nz/products/hygiene/careport/careport-freestanding-shower-cubicle/	Careport 120 Plus	All NZ	No	Yes	Yes	Shower only Internal size – 0.954mW x 1.156mD	<ul style="list-style-type: none"> • Freestanding portable shower unit • Addresses shorter-term needs • Toilet and/or basin can be added • Same specs as for Careport Extra